

# Infrared Analysis Request Form

Date:

Name:

Sample name or number:

Sample preparation:

pellet

mull

KBr plates

ATR

Scanwidth:

m.p. / b.p.

Special instructions:

## Instructions

1. Use a separate request form for each sample submitted.
2. A minimum of five milligrams of sample should be submitted.
3. Samples should be clearly labeled with your name and sample name or number.
4. Submit samples in screw - or snap-cap vials; no round bottom flasks will be accepted.
5. Place sample request form on the counter along with your sample (if submitting electronically, please place sample on the counter). The results can be picked up from the same place.
6. Radioactive samples will not be accepted.

Precautions:

explosive

volatile

hygroscopic

sensitive to \_\_\_\_\_

toxic, etc.